


PODIATRY AND MEDICAL CATALOG 2014



Equipped for Life™ with ALLPRO Imaging

Dear Valued Customer,

Since 1985, ALLPRO Imaging has provided the healthcare, security and non-destructive testing industries first-rate digital imaging products. Our mission continues to be strong: manufacture high quality imaging devices and software that simplify radiography procedures while providing excellent imagery.

As we progress into 2014, we promise to stand by our mission and bring innovation and variety to all of our customers. We are dedicated to each of our markets and will continue to strive for growth both domestically and internationally. ALLPRO Imaging has always been an innovative thinker, preparing for the changes and regulations occurring in our industries.

Thank you to our entire dealer network, partners, healthcare practitioners and industrial professionals for we could not continue to have success without your dedication and support. Our focus for the future will ensure that each of our customers are "Equipped for Life" as we continue to provide top-quality products and equipment.

Thank you again for all your support.

Regards,

A handwritten signature in black ink, appearing to read "David Chu". The signature is fluid and cursive, with the first name being more prominent than the last.

David Chu
Business Unit Manager



Podiatry



Digital Imaging

ScanX® 12 SE Packages	4
ScanX 12 SE Accessories	5
ScanX 12 Extended Warranties	5
ScanX Fit® Packages	6
ScanX Fit Accessories & Extended Warranties	7

General Medical



Digital Imaging

ScanX Fit® Packages	10
ScanX Fit Extended Warranties	11
ScanX Fit Accessories	13

Warranty Information



“The ScanX® 12 SE installation and training was quick and easy. My staff and I were up and running with this system in no time and we don’t miss our film processor at all!”

– Dr. Gregory Davies, DPM
Syosset, NY



PODIATRY

Digital Imaging

Accessories

Extended Warranties



scan-x

ScanX® 12 SE Digital Imaging Packages at a glance:

- Features high resolution digital images in less than 30 seconds.
- All plate sizes through 10" x 12"
- Hassle-free onsite installation and training, plus ongoing support during warranty period.
- Easy transition from film-based systems, minimal training needed by staff.
- Get up and running immediately; use your existing X-ray.
- Includes imaging software, computer workstation or laptop, imaging plates and accessories.



Packages Include:

ScanX 12 SE system with choice of laptop or desktop computer, with TigerView Software, 3 10" x 12" Imaging Plates, 3 10" x 12" Plate Protectors, 3 10" x 12" Cassettes. On-site installation. 1-Year Swap Warranty.

TigerView
Acquire | Store | View | Share

Part No. Description

G7660	ScanX 12 SE with Computer Workstation and TigerView Software (On-site Installation and Training)
G7670	ScanX 12 SE with Computer Laptop and TigerView Software (On-site Installation and Training)

ScanX 12 SE Specifications

Dimensions (WxDxH):

14"x 14"x 28" / 36cm x 36cm x 71cm

Weight:

55lbs. / 25Kgs.

Electrical:

60Hz, 115V

ScanX® 12 SE Package Accessories


Part No.	Description
B2010	Cleaning Sheets, Package of 12
B2020	Cleaning Sheets, Package of 25
B8910	PSP Cleaning Wipes, Box of 50
F3808	Dust Cover
73578-8M	Imaging Plate 8" x 10"
73578-10M	Imaging Plate 10" x 12"
74190-8	Cassette 8" x 10"
74190-10	Cassette 10" x 12"
B2160-8M	Plate Protector 8" x 10"
B2160-10	Plate Protector 10" x 12"



PSP Cleaning Wipes (PN B8910)

Extended Warranties

Part No.	Description
D5901-P	1-Year Extended Warranty for ScanX12 SE (Scanner Only)



Cassette 8" x 10" (PN 74190-8)
Cassette 10" x 12" (PN 74190-10)





ScanX Fit® Digital Imaging Packages at a glance:

- Includes user-friendly Imaging Software.
- Includes computer workstation or laptop, imaging plate and accessories.
- Hassle-free on-site installation and training by our experts, plus ongoing support while under warranty.
- Images stored in native DICOM 3.0 format; integration to EMR and Practice Management Systems with DICOM Modality Worklist.
- Get up and running immediately; use your existing x-ray system.
- Additional DICOM viewers functionality, and web server options, for access to images from any computer.
- Compact and portable system.



Packages Include:

ScanX Fit with choice of laptop or desktop computer, with TigerView Software, 3 10" x 12" Imaging Plates, 3 10" x 12" Plate Protectors, 3 10" x 12" Cassettes. On-site installation. 1-Year Swap Warranty.



Part No.	Description
G7620	ScanX Fit with Computer Workstation and TigerView Software (On-site Installation and Training)
G7640	ScanX Fit with Computer Laptop and Tigerview Software (On-site Installation and Training)

TigerView
Acquire | Store | View | Share

ScanX Fit Specifications

Dimensions (WxDxH):

15 1/2" x 18" x 14" / 39cm x 46cm x 36cm

Weight:

43lbs. / 19.5Kgs. (without battery)

45lbs. / 20.4Kgs. (with battery)

Electrical:

110-240 VAC, 50/60 Hz

ScanX Fit® Package Accessories

Part No.	Description
D5990	Hard Case
D5995	Soft Carrying Case
D5262	Dust Cover
D5940	10 Pack of Brushes
B2010	Cleaning Sheets, Package of 12
B2020	Cleaning Sheets, Package of 25
B8910	PSP Cleaning Wipes, Box of 50
73578-8M	Imaging Plate 8" x 10"
73578-10M	Imaging Plate 10" x 12"
73578-14M	Imaging Plate 14" x 17"
74190-8	Cassette 8" x 10"
74190-10	Cassette 10" x 12"
74190-14	Cassette 14" x 17"
B2160-8M	Plate Protector 8" x 10"
B2160-10	Plate Protector 10" x 12"
B2160-14	Plate Protector 14" x 17"



Plate Protectors 8" x 10" and 14" x 17" (PN B2160-8M and PN B2160-14M)

Extended Warranties

Part No.	Description
D5901-F	1-Year Extended Warranty for ScanX Fit (Scanner Only)



Soft Carrying Case (PN D5995)



“We chose the ScanX Fit® based on price, image quality and functionality for our practice. The training was excellent and our transition from film to digital was extremely smooth.”

Weaverville Family Medicine Associates
Weaverville, NC
Thomas Crook, Hospital Administrator
Ken Whitesides, Senior X-ray Technician


GENERAL MEDICAL



Digital Imaging

Accessories

Extended Warranties



scan-x fit

ScanX Fit® Digital Imaging Packages at a glance:

- Portable, compact system.
- Includes user-friendly TigerView Imaging Software.
- Customizable settings enable optimum image quality for your specific application.
- Includes computer workstation, imaging plates and accessories.
- Hassle-free on site installation and training by our experts, plus ongoing support while under warranty.
- Images stored in native DICOM 3.0 format; integration to EMR and Practice Management Systems with DICOM Modality Worklist.

Packages Include:

ScanX Fit with choice of laptop or desktop computer , TigerView DICOM Software with Image Acquisition and 2 DICOM Viewer licenses, 2 14" x 17" Imaging Plates, 2 14" x 17" Cassettes, 2 14" x 17" Plate Protectors, On-site installation. 1-Year Standard Swap Warranty.



Part No.	Description
G7620-M	ScanX Fit with Computer Workstation and Tigerview Software (On-site Installation and Training)
G7640-M	ScanX Fit with Computer Laptop and Tigerview Software (On-site Installation and Training)

TigerView

Acquire | Store | View | Share

Specifications

Dimensions (WxLxH):

15 ½" x 18" x 14" / 39cm x 46cm x 36cm

Weight:

43lbs. / 19.5Kgs. (without battery)

45lbs. / 20.4Kgs. (with battery)

Electrical:

110-240 VAC, 50/60 Hz

Extended Warranties

Part No.	Description
D5901-F	1-Year Extended Warranty for ScanX Fit (Scanner Only)





ScanX Fit® Package Accessories

Part No.	Description
D5990	Hard Case
D5995	Soft Carrying Case
D5940	10 Pack of Brushes
B2010	Cleaning Sheets, Package of 12
B2020	Cleaning Sheets, Package of 25
B8910	PSP Cleaning Wipes, Box of 50
73578-8M	Imaging Plate 8" x 10"
73578-10M	Imaging Plate 10" x 12"
73578-14M	Imaging Plate 14" x 17"
74190-8	Cassette 8" x 10"
74190-10	Cassette 10" x 12"
74190-14	Cassette 14" x 17"
B2160-8M	Plate Protector 8" x 10"
B2160-10	Plate Protector 10" x 12"
B2160-14	Plate Protector 14" x 17"




Imaging Plates 8" x 10" and 14" x 17"
(PN 73578-8M and PN 73578-14M)



PSP Cleaning Wipes (PN B8910)



Soft Carrying Case (PN D5995)



Imaging Cassettes 8" x 10" and 14" x 17"
(PN 74190-10 and PN 74190-14M)



Hard Case (PN D5990)




ALL I MA



TigerView
Acquire | Store | View | Share

PRO AGING





Standard Warranty and Returns Policy

All parts and component returns and replacement of equipment under warranty require a Return Materials Authorization (RMA) issued by Tech Service.

Warranty credit will not be applied to product returns that exhibit damage due to shipping, misuse, careless handling or repairs by unauthorized personnel. Credit, or partial credit, will not be issued until product/parts have been received and assessed. Warranty is void if product is installed or serviced by anyone other than authorized ALLPRO Imaging dealer service personnel.

Dealer labor, shipping & handling charges are not covered by this warranty.

ALLPRO Imaging will institute a 25% restocking fee. Spare parts and accessories older than 12 months from invoice date and/or products with a tamper-resistant seal (including electrical products, PSP, barrier envelopes, sheaths, etc.) that have been opened, cannot be returned.

ALLPRO Imaging equipment is warranted to be free from defects in material and workmanship from the date of installation for a period of 12 months.



Out-of-Warranty Programs

For information on our S.U.R.E. Out-of-Warranty, Swap*, Upgrade, Trade-In and Service options, please call 888-862-4050 or visit www.allproimaging.com/healthcare to download S.U.R.E. Program literature.

*The Swap Warranty is valid in the US and Canada only.

There is a minimum order requirement of \$75.

Contact Information

For general and immediate questions, please call:

Toll Free in the US and Canada: 1-800-AIR TECH (247-8324)
Monday through Friday, 8:00 AM – 8:00 PM (EST)

Disclaimer

1. Items returned without an RMA, or included with other products for which an RMA has been issued, may be returned to the customer at the discretion of ALLPRO Imaging.
2. Advanced replacement will be provided and invoiced at full price.
3. For Shipped Items – ALL freight costs are paid for by dealers, unless part of an approved promotion/program. All shipping damage costs should be reviewed and addressed with carrier or customer/dealer as needed.
4. \$75 minimum order size for all non-warranty orders.



ALLPRO Imaging is a leading manufacturer of imaging systems dedicated to the medical, veterinary, security and industrial markets for over 25 years. ALLPRO maintains a 200,000 square-foot manufacturing facility in Melville, New York, and a network of authorized dealers.

To locate your authorized dealer call **888.862.4050**.

1295 Walt Whitman Road | Melville, NY 11747-3062 | Phone 888.862.4050 | Fax 888.247.8481

www.allproimaging.com



**ALLPRO
IMAGING**
equipped for life™

PN 9922-134 Rev. J