



2006

Parts Price List

Suggested Retail Prices
Effective January 1, 2006

This price list supersedes all other lists.
All prices are subject to change without notice.
Shipments are F.O.B., Hicksville, New York.

TABLE OF CONTENTS

RETURNS	4
REBUILT PARTS EXCHANGE POLICY	4
PARTS WARRANTY POLICY	4
CONTACT INFORMATION	5
RETURN AUTHORIZATION POLICY (RA)	5
ALLPRO 100 PLUS	6
ALLPRO 100 PLUS	7
ALLPRO 100 PLUS	8
ALLPRO 100 PLUS EXCHANGE PARTS	8
AP2010 BASE ASSEMBLY	9
AP2010 DRIVE BLOCK ASSEMBLY	10
AP2010 COVER ASSEMBLY	10
ALLPRO 2000M	11
ALLPRO 2000M COVER ASSEMBLY	12
ALLPRO 2000M BASE ASSEMBLY	12
ALLPRO 2000M DRIVE BLOCK ASSEMBLY	13
AP2010/AP2000M DEVELOPER RACK	13
AP2010 / AP2000M FIXER RACK	14
AP2010 / AP2000M WASH/DRY ASSEMBLY	14
ELECTRONIC ASSEMBLY, AUTOMATIC PROCESS	15
ELECTRONIC ASSEMBLY, MANUAL PROCESS	15
10 INCH RACKS	15
FILM TRANSPORT (8 INCH)	15
ACCENT CCD DIRECT DIGITAL INTRAORAL SENSOR	16
AP SCANX	17

PARTS WARRANTY POLICY

Should any part prove defective within three months from the date of installation, it will be repaired or replaced, without charge, upon receipt at our factory in Hicksville, New York, transportation charges prepaid.

Our obligation under this warranty is limited to putting in proper operating condition any part or parts returned to us through the authorized dealer from whom originally purchased and which our examination shall disclose to have been defective.

This warranty does not apply to damage due to misuse or careless handling. ALLPRO Imaging is not liable for indirect or consequential damage or loss of any nature in connection with parts sold by them. No charge for labor or expense is required to repair defective goods or damage occasioned by them will be allowed.

This warranty is in lieu of all other warranties, expressed or implied, and no representative or person is authorized to assume for us any other liability in connection with the sale of parts.

REBUILT PARTS EXCHANGE POLICY

Certain parts are available on an "Exchange Basis." Rebuilt "Exchange" parts may be ordered in either of two ways:

1. Advance Replacement

- A) Order the rebuilt part needed.
- B) We ship and invoice at the "Outright Sale Price."
- C) Return the old part (freight prepaid) to factory.
- D) We issue the credit for difference between "Outright Sale Price" and "Exchange Price."

2. Direct Exchange

- A) Return the old part (Freight Prepaid) to the factory.
- B) Request that we replace it with a rebuilt part.
- C) We ship the rebuilt part (Freight Collect).
- D) We invoice at the "Exchange Price."

RETURNS

If ordered in error, goods may be returned for a credit/refund subject to the following conditions:

- 1. All returned goods require a Return Authorization (RA) number.
- 2. Goods must be returned within 30 days of shipment.
- 3. Goods must be returned undamaged and unused in original containers.
- 4. Transportation charges to and from the purchaser must be paid by the purchaser.
- 5. A 15% restocking charge is applied toward returned goods.

RETURN AUTHORIZATION POLICY (RA)

All part and component returns and replacement of equipment under warranty require a Return Authorization (RA) number.

Please follow this procedure:

1. Call Technical Service or Customer Service. The service representative will issue you a RA number after you provide him with the following:
 - Model Number
 - Serial Number of the piece of equipment that needs repair
 - Part number and serial number (when applicable) of the components being returned
2. Record the RA number on the paperwork you are enclosing and on the outside of the shipping container before returning it to ALLPRO Imaging.
3. Only these managers can issue a RA number for equipment:
 - Technical Service Manager
 - Customer Service Manager
 - National Sales Manager

CONTACT INFORMATION

Our Customer Relations Department is available:
Monday - Friday, 8:30 AM - 5 PM E.S.T.;

Phone: 1-800-247-8324 or 1-516-433-7676

Fax: 1-888-247-8481 or 1-516-433-7684

ALLPRO 100 PLUS

PART NO.	DESCRIPTION	PRICE
47465	ACTUATOR,FILM SENSOR	\$3.15
49816	BEARING, RACK DRIVE SHAFT	\$1.30
49241	BELT,MAIN DRIVE,133-TOOTH(100 SEC.OPER.)(FRONT MOTOR REMOVAL)	\$22.80
47124	BRACKET, WASH TUBE SUPPORT	\$14.10
47123	BRACKET,FEEDTRAY HOLD DOWN	\$7.05
47318	BRACKET,HIGH LEVEL SENSOR	\$2.00
49124	BRACKET,SIDE PANEL	\$3.60
49811	BUSHING, ROLLER , SMALL	\$3.15
49831	BUSHING, ROLLER, LARGE	\$3.15
49161	CATCH,BULLET TYPE,SIDE PANEL	\$1.60
49511	CIRCUIT BREAKER, 115V	\$39.70
30800M	COLLAR	\$6.90
47182	COUPLING, LEFT	\$7.90
47181	COUPLING,RIGHT	\$7.90
49727	COVER PLATE, FRONT	\$43.15
47185	COVER,FILM SENSOR	\$7.95
47184	COVER,INTERLOCK SWITCH	\$11.75
49126	COVER,MAIN CONTROL BOARD, TRANSPARENT	\$75.05
49713	CROSS TIE, DRYER RACK	\$9.80
49827	CROSS TIE, WET TRANSPORT	\$9.80
47712	DEVELOPER CHEMISTRY FLOATING LID	\$15.10
47053	DEVELOPER REPLENISHMENT CONTAINER ASSY	\$61.85
47810	DEVELOPER TRANSPORT ASSEMBLY	\$1,175.00
47028	DIAGNOSTIC DOOR ASSEMBLY	\$10.75
49461	DRIVE MOTOR ASSY,115V (MOUNTING PLATE,DRIVE MOTOR, PULLEY & FAN)	\$253.90
49467	DRIVE MOTOR ASSY,230V (MOUNTING PLATE,DRIVE MOTOR,PULLEY & FAN)	\$273.30
49301	DRIVE MOTOR KIT, 115V	\$200.15
49312	DRIVE MOTOR KIT, 230V	\$223.20
49737	DRYER AIR TUBE, LONG, RIGHT SIDE	\$9.25
49738	DRYER AIR TUBE, SHORT, LEFT SIDE	\$9.25
49846	DRYER BLOWER, 115V	\$180.90
47718	DRYER SIDE PLATE, LEFT	\$139.05
47717	DRYER SIDE PLATE, RIGHT	\$139.05
47750	DRYER TRANSPORT ASSEMBLY, 115V	\$2,225.00
47765	DRYER TRANSPORT ASSEMBLY, 230V	\$2,225.00
47737	DRYER, REAR SUPPORT ASSEMBLY	\$15.20
47586	ELECTRONIC ENCLOSURE ASSEMBLY, 230V	\$1,740.70
47520	ELECTRONICS ENCLOSURE ASSEMBLY, 115V	\$1,735.55
47190	FEED TRAY ASSEMBLY	\$288.40
47798	FILM SENSOR ACTUATOR ASSEMBLY	\$14.95
49195	FILM SENSOR WIRE HARNESS	\$9.50
47054	FIXER REPLENISHMENT CONTAINER ASSY	\$61.85
47820	FIXER TRANSPORT ASSEMBLY	\$795.00
47188	GASKET, TOP COVER	\$4.55
49802	GEAR, IDLER, 1 1/4" DIA.	\$6.55
49804	GEAR, IDLER, 2.0" DIA.	\$7.20
49818	GEAR, MAIN DRIVE GEAR, MAIN DRIVE SHAFT SIDE (GRAY)	\$9.50
49803	GEAR, MAIN DRIVE, GEAR SIDE, 1 1/2" DIA.	\$6.60
49801	GEAR, ROLLER, 3/4" DIA.	\$1.55
49805	GEAR, SPUR, .75P.D., 3/8" BORE	\$3.90
47902	GEAR,20 TOOTH	\$10.15
49711	GROMMET,5/16".I.D.	\$1.65
49616	GUARD, LAMP, SNAP-ON	\$25.75
47768	GUIDE, CROSSOVER, DRYER INLET, FILM TRANSPORT	\$28.35
47869	GUIDE, CROSSOVER, WASH OUTLET, FILM TRANSPORT	\$26.30
49814	GUIDE, CURVED	\$29.50

ALLPRO 100 PLUS CONTINUED

PART NO.	DESCRIPTION	PRICE
49824	GUIDE, CURVED, ARCHED	\$27.60
49815	GUIDE, FLAT	\$33.80
47868	GUIDE,CROSSOVER,IN/OUT,FILM TRANSPORT	\$27.30
47569	HARNESS, MAIN ASSEMBLY	\$189.15
49408	HEAT SINK ASSEMBLY, 115V, INCLUDES TRIACS	\$101.40
49336	HEATER, DEVELOPER, 115V	\$182.10
49575	HEATER, DEVELOPER, 230V	\$182.10
47519	INTERLOCK SWITCH	\$18.55
49142	J-TUBE,REPLENISHER & WATER TANK INLET	\$6.10
47143	KIT:DRAIN VALVE WASH/WEEP ASSEMBLY	\$57.80
47980	KIT;FILM SENSOR PCB REPLACEMENT	\$34.00
49789	LAMP GUARD	\$66.75
49714	LAMP HOLDER	\$2.15
49751	LAMP, INFRARED, FRONT, 115V DRYER	\$133.50
49253	LAMP, INFRARED, FRONT, 230V	\$131.45
49734	LAMP, INFRARED, MIDDLE, 115V - DRYER	\$133.50
49252	LAMP, INFRARED, MIDDLE, 230V	\$131.45
49752	LAMP, INFRARED, REAR, 115V DRYER	\$133.50
49254	LAMP, INFRARED, REAR, 230V	\$131.45
43106	LEVELING FOOT	\$1.95
49206	LOCKING TAB,SIDE PANEL	\$5.15
47103	MANUAL & BINDER, 100 PLUS - 115V	\$71.55
47104	MANUAL & BINDER, 100 PLUS - 230V	\$71.55
47725	PANEL, FRONT BAFFLE	\$46.25
47728	PANEL,DRYER BACK	\$113.30
30799	PIN,COTTER 1 1/2"	\$0.10
30798	PIN,COTTER 1"	\$0.05
117488	PLUG, FLOOD REPLENISHMENT (MAIN CONTROL BOARD)	\$5.00
49411	POWER CORD ASSEMBLY	\$32.35
49417	PULLEY, MAIN DRIVE MOTOR, 28-TOOTH (100 SEC.OPERATION)	\$22.40
49418	PULLEY, MAIN DRIVE SHAFT, 24-TOOTH	\$21.80
47314	RECIRCULATOR PUMP, DUAL HEAD, 115V	\$406.85
47304	RECIRCULATOR PUMP, DUAL HEAD, 230V	\$406.85
47774	REFLECTOR, BOTTOM	\$38.90
47721	REFLECTOR, SIDE	\$37.00
47722	REFLECTOR, TOP	\$28.00
49880	REPLENISHER PUMP ASSEMBLY, 115V (IWAKI PUMP)	\$173.35
30840	RETAINING CLIP,E-TYPE	\$1.30
49872	ROLLER, PLASTIC, GOLDEN COLOR	\$45.45
49828	ROLLER, RUBBER	\$72.45
49835	ROLLER, RUBBER, 20 - 30 MICROINCH	\$99.80
49307	SADDLE BEARING FOR TRANSPORTS & MAIN DRIVE SHAFT	\$3.70
30405	SCREW, SIDE PANEL LOCK	\$0.20
49817	SHAFT DRIVE 3/8 DIA.	\$20.05
49741	SHAFT DRIVE,DRYER	\$22.05
47230	SIDE PANEL ASSEMBLY,LEFT	\$188.70
47220	SIDE PANEL ASSEMBLY,RIGHT	\$188.70
47836	SIDE PLATE, DEVELOPER TRANSPORT	\$19.00
47837	SIDE PLATE, FIXER/WASH TRANSPORT	\$15.90
47848	SLEEVE,LEVELING FOOT	\$0.05
49808	SNAP RING, LARGE	\$0.20
49807	SNAP RING, SMALL	\$0.20
47450	SOLENOID VALVE ASSEMBLY, 115V, WATER INLET	\$86.50
47240	SOLENOID VALVE ASSEMBLY, 230V, WATER INLET	\$100.25
47199	SPLASH GUARD	\$99.80
47163	SPLASH GUARD ASSEMBLY	\$204.60

ALLPRO 100 PLUS CONTINUED

PART NO.	DESCRIPTION	PRICE
49791	SPRING, LAMP GUARD RETAINER	\$1.35
49822	STUD, IDLER	\$5.55
47302	TANK, DRILLED	\$278.10
49823	TENSION BAND, COLOR CODED	\$1.55
47335	THERMISTOR ASSEMBLY, DEVELOPER	\$111.00
49736	THERMOSTAT (OVER-TEMPERATURE PROTECTION)	\$44.65
47110	TOP COVER ASSEMBLY	\$437.75
115578	TRIAC,600V,25 A	\$51.90
47614	TUBE ASSY,DEVELOPER REPLENISHER	\$15.70
47615	TUBE ASSY,FIXER REPLENISHER	\$15.70
47141	VALVE, DEVELOPER DRAIN ASSEMBLY	\$57.80
47142	VALVE, FIXER DRAIN ASSEMBLY	\$57.80
49102	WALL CHART	\$4.70
47830	WASH TRANSPORT ASSEMBLY	\$860.00
49166	WATER INLET HOSE, RIGHT ANGLE	\$29.30
49730	HARNES DRYER ASSEMBLY	\$111.25
49120	WORM SHAFT ASSEMBLY, COMPLETE	\$191.00
49123	WORM WITH BUSHING	\$14.00
49195	HARNES SHUTTER CONNECTOR ASSEMBLY	\$9.50

ALLPRO 100 PLUS EXCHANGE PARTS

PART NO.	DESCRIPTION	OUTRIGHT SALE PRICE	EXCHANGE SALE PRICE
47750	DRYER ASS'Y, 115V	\$2,225.00	\$1,231.00
47765	DRYER ASSY, 230V	\$2,225.00	\$1,231.00
47810	DEVELOPER RACK ASS'Y	\$1,175.00	\$654.00
47500	MAIN CONTROL BOARD, 115V	\$1,125.00	\$494.00
47582	MAIN CONTROL BOARD, 230V	\$1,125.00	\$494.00
47820	FIXER RACK ASS'Y	\$795.00	\$453.00
47830	WASH RACK ASS'Y	\$860.00	\$479.00

AP2010 BASE ASSEMBLY

PART NO.	DESCRIPTION	PRICE
45201	BASE MOLDED W/HARDWARE INSERTS & DRAIN ELBOWS	\$587.10
45440	BASE WIRING HARNESS, COMPLETE	\$132.75
43232	BEARING CERAMIC	\$23.40
43885	BOTTOM PLATE ASSEMBLY	\$42.45
45262	CIRCULATOR DRAW ASSEMBLY	\$272.95
45691	CIRCULATOR DRAWER REPLACEMENT KIT (INCLUDES MAGNETS)	\$334.75
49963	CLAMP, TANK DRAIN	\$1.25
45445	CONNECTOR PLATE, ETCHED	\$15.65
43384	DEVELOPER CIRCULATOR COVER - BLACK	\$6.60
43357	DRAIN HOSE STRAP	\$4.55
44361	DRAIN HOSE, DEVELOPER (BLACK)	\$18.20
44362	DRAIN HOSE, FIXER (RED)	\$18.20
44363	DRAIN HOSE, WATER (WHITE)	\$18.20
44206	DRAIN MANIFOLD, MOLDED	\$21.95
43286	DRAIN TUBE ASSY - BLACK	\$13.25
43287	DRAIN TUBE ASSY - RED	\$13.25
43270	DRAIN TUBE ASSY - TAN	\$11.90
43271	DRAIN TUBE WASHER	\$1.35
45704	DRIVE BLOCK ASSEMBLY AP2010	\$973.35
43259	ELECTRICAL COMPARTMENT COVER	\$53.25
45975	EMI/RFI FILTER REPLACEMENT, IEC POWER ENTRY	\$34.10
45436	EXTERIOR IEC MOUNTING PLATE	\$1.10
45567	FILM INLET MODIFIED	\$20.75
43385	FIXER CIRCULATOR COVER - RED	\$6.65
90219	FUSE HOLDER	\$7.20
45713	FUSE, DELAYED ACTION 10 AMP	\$1.35
43258	GEAR COVER, MOLDED	\$40.80
43859	HEATER PAD HOLD DOWN	\$6.65
45437	INTERIOR IEC MOUNTING PLATE	\$2.80
45100	LEVELING BASE ASSEMBLY	\$154.60
43106	LEVELING FEET	\$1.95
45920	LINE CORD REPLACEMENT, 115V PLUG IEC SOCKET	\$12.20
45950	PAD ASSY REPLACEMENT	\$201.20
49793	POPPET VAVE KIT (INCLUDES 2 VALVES & 2 O-RINGS)	\$50.00
44265	REPLENISHER ASSEMBLY	\$357.05
45935	REPLENISHER ASSY W/CIRC.DRAWER ASSY	\$629.35
44044	REPLENISHER DEVELOPER INLET TUBING (BLK)	\$5.00
44042	REPLENISHER DEVELOPER OUTLET TUBING (BLK)	\$3.00
44041	REPLENISHER FIXER OUTLET TUBING (RED)	\$1.90
44043	REPLENISHER FIXER OUTLET TUBING (RED)	\$5.00
44039	REPLENISHER TUBING CLAMP	\$1.90
30569	SCREW - DRIVE BLOCK MOUNT	\$0.05
30539	SCREW - GEAR COVER	\$0.05
30568	SCREW - GEAR COVER	\$0.10
30550	SCREW, BOTTOM PLATE	\$0.20
45980	SOLENOID ASSEMBLY	\$74.60
45323	STANDOFF, INTERLOCK ACTUATOR	\$5.10
41286	TEFLON WASHER	\$0.20
45970	TEMP SENSOR	\$43.80
30552	THUMB SCREW, CIRCULATOR COVER	\$0.25
44210	TUBING FIXER DRAIN	\$2.50
44211	TUBING WATER DRAIN	\$1.75
44209	TUBING, DEVELOPER DRAIN	\$6.75
43189	WASH TANK DRAIN TUBE PLUG KIT	\$2.20
45221	WATER INLET J TUBE	\$3.35
45450	WATER LEVEL SENSOR	\$40.45

AP2010 DRIVE BLOCK ASSEMBLY

PART NO.	DESCRIPTION	PRICE
43244	BEARING, SHAFT	\$1.70
44713	CAPACITOR;DRIVE MOTOR	\$19.40
45704	DRIVE BLOCK ASSY	\$973.35
43247	DRIVE CHAIN	\$30.55
43243	DRIVE GEAR (12 TOOTH)	\$4.05
43253	DRIVE MOTOR MOUNTING PLATE	\$3.65
43246	DRIVE SHAFT CLAMP - RED	\$1.70
43223	DRIVE SHAFT CLAMP- BLACK	\$1.70
43242	MAIN DRIVE SHAFT	\$72.50
43243-1	MOTOR DRIVE GEAR 16 TOOTH	\$4.05
44715	PCB ASSY; BASE, INTERCONNECT	\$69.55
43254	RACK LATCH - BLACK	\$2.80
41127	SNAP RING	\$1.75
45970	THERMISTOR ASSY	\$43.80
45450	WATER LEVEL SENSOR	\$40.45

AP2010 COVER ASSEMBLY

PART NO.	DESCRIPTION	PRICE
43605	AIR HEATER ONLY, 600W	\$109.10
43672	AIR HEATER REPLACEMENT KIT	\$113.30
43797-1	CABLE ASSY;MEDICAL FAN PANEL	\$20.05
45350	COVER INSERT, SILKSCREENED	\$56.55
45910	COVER INTERLOCK SWITCH KIT	\$40.85
43868	DRYER ASSEMBLY, COMPLETE	\$535.60
45715	ELECTRONICS MODULE	\$2,270.00
44482	FILM SENSOR ASSEMBLY	\$32.85
43587	HINGE BLOCK	\$5.95
45360	MOLDED COVER W/HARDWARE INSERTS	\$367.70
43584	SHUTTER ASSEMBLY	\$26.25
43659	SHUTTER GUARD	\$25.00

ALLPRO 2000M

PART NO.	DESCRIPTION	PRICE
49794	BELLOWS, REPLENISHER PUMP FOR P/N 43045	\$118.00
43754	BELLOWS, REPLENISHER PUMP, FOR P/N 40703-1	\$34.60
43218	CLAMP	\$3.10
40153	CUFF ASSY FOR DAYLIGHT LOADER	\$186.60
43513	DRIVE GEAR (14-TOOTH)	\$4.00
43281	DRIVE MOTOR REPLACEMENT KIT	\$488.00
43403	ELECTRONICS MODULE (MANUAL PROCESS)	\$2,270.00
40719	FUSE, 2 AMP (PACKAGE OF 5)	\$5.95
43503	GEAR (28-TOOTH)	\$1.45
43504	GEAR (36-TOOTH)	\$1.45
43520	GEAR (60-TOOTH)	\$3.70
43736	GEAR REPLACEMENT KIT; LG GEARS	\$10.10
43290	HEATER BAR ASSEMBLY	\$405.10
43694	IDLER GEAR REPLACEMENT KIT	\$9.20
40782	INSERT, STAINLESS STEEL	\$2.40
43869	KEY FOR DRIVE MOTOR BOLTS	\$8.55
43942	KIT, DRIVE MOTOR PCB & HARNESS ASSY	\$73.50
40715	NUT, NYLON - RED	\$1.25
43490	PC BOARD, AUTOMATIC REPLENISHER	\$426.40
43498	PC BOARD, AUTOMATIC REPLENISHER (EXT. ADJ)	\$422.30
43420	PC BOARD, BASE (D/F HEATER BAR)	\$121.55
43771	POWER CABLE, REPLENISHER	\$117.85
43090	REMOTE REPLENISHER CONTROL	\$98.30
43720	REPLENISHER BACK PANEL ASSEMBLY	\$959.45
43045	REPLENISHMENT PUMP AND MOTOR ASSY	\$519.90
41292	ROLLER GEAR, 3/16" THICK (14-TOOTH)	\$1.35
43842	ROLLER, PLASTIC (10 INCH)	\$26.80
43518	ROLLER, PLASTIC (8 INCH)	\$26.80
43843	ROLLER, RUBBER (10 INCH)	\$58.10
43519	ROLLER, RUBBER (8 INCH)	\$51.90
43509	SHAFT GEAR (24-TOOTH)	\$3.70
43847	SPACER BAR (10 INCH)	\$15.65
41285	SPACER BAR (8 INCH)	\$6.05
43517	SPACER RETAINER	\$2.40
43524	SPACER WASHER	\$3.05
40727	TUBING, BLACK (8 FEET)	\$5.00
40726	TUBING, RED (8 FEET)	\$5.00
49793	VALVE, REPLENISHER PUMP (2 EA) FOR P/N 43045	\$50.00
43752	VALVE, REPLENISHER PUMP, FOR P/N 40703-1	\$45.40
43543	VARISTOR ASSY	\$7.40
41320	WASHER, BLACK	\$1.30
41821	WASHER, FLAT (BLACK)	\$0.10
41284	WASHER, FLAT (WHITE)	\$0.20
41248	WASHER, RUBBER	\$1.40

ALLPRO 2000M COVER ASSEMBLY

PART NO.	DESCRIPTION	PRICE
43805M	COVER ASSEMBLY WITH ELECTRONICS	\$3,288.80
43803	COVER ASSEMBLY WITH ELECTRONICS (MANUAL PROCESS)	\$3,288.80
43530	DRYER FAN	\$314.50
43715M	ELECTRONICS ASSEMBLY	\$2,270.00
43868	FAN AND DRYER ASSEMBLY, 600W	\$535.60
43605	HEATER, DRYER, 600W	\$109.10
43587	HINGE BLOCK ASSEMBLY	\$5.95
43495	PC BOARD AND CABLE ASSY, SHUTTER	\$175.90
43459	POWER SWITCH	\$14.15
43477	PROCESS SWITCH (AUTOMATIC PROCESS)	\$14.15
43469	PROCESS SWITCH (MANUAL PROCESS)	\$13.75
43115	PROP STICK	\$9.20
43584	SHUTTER ASSEMBLY	\$26.25
43659	SHUTTER GUARD ASSEMBLY, MEDICAL	\$25.00

ALLPRO 2000M BASE ASSEMBLY

PART NO.	DESCRIPTION	PRICE
43371	90 DEG. WATER INLET TUBE	\$24.00
43802	BASE ASSEMBLY (D/F HEATER BAR)	\$3,218.75
43804	BASE ASSEMBLY COMPLETE (D/F HEATER PAD)	\$3,218.75
43232	BEARING, CERAMIC	\$23.40
43267	CIRCULATOR DRIVE MAGNET	\$35.30
43233	CIRCULATOR MAGNET	\$27.60
43695	CIRCULATOR MOTOR, FAN AND MAGNET ASSEMBLY	\$88.85
41211	CLAMP, DRAIN HOSE	\$3.70
43236	CLAMP, HOSE	\$2.60
43384	DEVELOPER CIRCULATOR COVER (BLACK)	\$6.60
41212	DRAIN HOSE	\$46.00
44206	DRAIN MANIFOLD ASSEMBLY; NEW STYLE	\$21.95
43286	DRAIN TUBE ASSEMBLY - BLACK	\$13.25
43287	DRAIN TUBE ASSEMBLY - RED	\$13.25
43270	DRAIN TUBE ASSEMBLY - TAN	\$11.90
43271	DRAIN TUBE WASHER	\$1.35
43240	DRIVE BLOCK ASSEMBLY	\$922.15
43385	FIXER CIRCULATOR COVER (RED)	\$6.65
43540	FREQUENCY FILTER AND VARISTOR ASSEMBLY	\$58.10
90219	FUSE HOLDER	\$7.20
41914	FUSE, SLO-BLOW, 8AMP (PACKAGE OF 5)	\$19.00
43258	GEAR COVER	\$40.80
43864	HARNESS, BASE	\$105.50
43863	HARNESS, BASE (D/F HEATER BAR)	\$96.40
43858	HEATER PAD	\$238.80
43859	HEATER PAD HOLD-DOWN	\$6.65
49165	HOSE ASSEMBLY, WATER (AMERICAN)	\$26.70
44361	HOSE, DEVELOPER DRAIN	\$18.20
44362	HOSE, FIXER DRAIN	\$18.20
44363	HOSE, WASH DRAIN	\$18.20
43435	LINE CORD	\$8.95
43206	MANIFOLD & HOSE ASSEMBLY	\$66.35
40806	NUT, NYLON (BLACK)	\$1.25
43362	O-RING - HEATER BAR	\$1.45
43485	PC BOARD, BASE (D/F HEATER PAD)	\$90.65

ALLPRO 2000M BASE ASSEMBLY CONTINUED

PART NO.	DESCRIPTION	PRICE
43699	RECOVERY DRAIN HOSE KIT	\$123.85
43225	SOLENOID VALVE ASSEMBLY	\$115.15
43214	TANK INLET DEVELOPER REPLENISHER TUBING	\$2.80
43216	TANK INLET FIXER REPLENISHER TUBING	\$1.90
43856	TEMPERATURE SENSOR ASSEMBLY	\$62.70
30552	THUMBSCREW, CIRCULATOR COVER	\$0.25
43866	TUBE, WATER INLET 22"	\$6.95
43447	WATER LEVEL SENSOR	\$40.45

ALLPRO 2000M DRIVE BLOCK ASSEMBLY

PART NO.	DESCRIPTION	PRICE
43244	BEARING, SHAFT	\$1.70
43250	BRACKET,PC BOARD DRIVE MOTOR	\$3.75
110335	CAPACITOR, DRIVE MOTOR	\$25.25
43247	DRIVE CHAIN	\$30.55
43243-2	DRIVE GEAR (10 TOOTH)	\$8.45
43243	DRIVE GEAR (12 TOOTH)	\$4.05
43243-1	DRIVE GEAR (16 TOOTH)	\$4.05
43248	DRIVE MOTOR ASSY (INCLUDES P/N 43178, 43281,41127,43942)	\$577.60
43253	DRIVE MOTOR MOUNT	\$3.65
43223	DRIVE SHAFT CLAMP - BLACK	\$1.70
43246	DRIVE SHAFT CLAMP - RED	\$1.70
43661	KEY FOR OPTICAL WHEEL SET SCREW	\$7.15
43242	MAIN DRIVE SHAFT	\$72.50
43178	OPTICAL SENSOR WHEEL REPLACEMENT KIT	\$37.65
43942	PCB & HARNESS ASSY, DRIVE MOTOR	\$73.50
43254	RACK LATCH - BLACK	\$2.80
30687	SET SCREW, OPTICAL SENSOR WHEEL	\$0.45
41127	SNAP RING	\$1.75

AP2010/AP2000M DEVELOPER RACK

PART NO.	DESCRIPTION	PRICE
43508	AXLE, SHAFT	\$4.70
43513	DRIVE GEAR (14 TOOTH)	\$4.00
43502	END PLATE	\$59.35
43503	GEAR (28 TOOTH)	\$1.45
43504	GEAR (36 TOOTH)	\$1.45
43736	GEAR REPLACEMENT KIT-(28,36,60 TEETH GEARS)	\$10.10
43846	GUIDE ROLLER	\$32.50
43516	HANDLE	\$7.90
41292	IDLER GEAR (14 TOOTH) 3/16" THICK	\$1.35
43694	IDLER GEAR REPLACEMENT KIT	\$9.20
43518	PLASTIC ROLLER	\$26.80
43843	RUBBER ROLLER	\$58.10
41284	RUBBER ROLLER SPACER, WHITE	\$0.20
30567	SCREW	\$0.10
30595	SCREW	\$0.15
43509	SHAFT GEAR (24 TOOTH)	\$3.70
43524	SPACE WASHER	\$3.05
41285	SPACER BAR	\$6.05
43517	SPACER RETAINER	\$2.40

AP2010 / AP2000M FIXER RACK

PART NO.	DESCRIPTION	PRICE
43508	AXLE SHAFT	\$4.70
43513	DRIVE GEAR (14 TOOTH)	\$4.00
43507	END PLATE	\$63.80
43503	GEAR (28 TOOTH)	\$1.45
43504	GEAR (36 TOOTH)	\$1.45
43736	GEAR REPLACEMENT KIT (28,36,60 TEETH GEARS)	\$10.10
43846	GUIDE ROLLER	\$32.50
43516	HANDLE	\$7.90
41292	IDLER GEAR (14 TOOTH) 3/16" THICK	\$1.35
43694	IDLER GEAR REPL KIT - 14 TOOTH	\$9.20
43842	PLASTIC ROLLER	\$26.80
41821	RUBBER ROLLER SPACER BLACK	\$0.10
41284	RUBBER ROLLER SPACER, WHITE	\$0.20
30567	SCREW	\$0.10
30595	SCREW	\$0.15
43509	SHAFT GEAR (24 TOOTH)	\$3.70
41285	SPACER BAR	\$6.05
43517	SPACER RETAINER	\$2.40
43524	SPACER WASHER	\$3.05

AP2010 / AP2000M WASH/DRY ASSEMBLY

PART NO.	DESCRIPTION	PRICE
43508	AXLE, SHAFT	\$4.70
43513	DRIVE GEAR (14 TOOTH) :	\$4.00
43503	GEAR (28 TOOTH)	\$1.45
43504	GEAR (36 TOOTH)	\$1.45
43520	GEAR (60 TOOTH)	\$3.70
43694	GEAR REPLACEMENT KIT (14 TOOTH GEAR)	\$9.20
43736	GEAR REPLACEMENT KIT (28,36,60 TEETH GEARS)	\$10.10
43846	GUIDE ROLLER	\$32.50
41292	IDLER GEAR (14 TOOTH) 3/16" THICK	\$1.35
41307	INNER AIR BAFFLE	\$11.40
43526	OUTER AIR BAFFLE	\$52.50
43842	PLASTIC ROLLER	\$26.80
43512	RIGHT END PLATE	\$105.70
43843	RUBBER ROLLER	\$58.10
41821	RUBBER ROLLER SPACER BLACK	\$0.10
41284	RUBBER ROLLER SPACER WHITE	\$0.20
30550	SCREW	\$0.20
30567	SCREW	\$0.10
43509	SHAFT GEAR (24 TOOTH)	\$3.70
41285	SPACER BAR	\$6.05
43517	SPACER RETAINER	\$2.40
43524	SPACER WASHER	\$3.05

ELECTRONIC ASSEMBLY, AUTOMATIC PROCESS

PART NO.	DESCRIPTION	OUTRIGHT SALE PRICE	EXCHANGE SALE PRICE
43626	FILM QUICK CT	\$2,270.00	\$942.00
43715M	ALL-PRO 2000M	\$2,270.00	\$942.00
45715	ALL-PRO 2010	\$2,270.00	\$942.00

ELECTRONIC ASSEMBLY, MANUAL PROCESS

PART NO.	DESCRIPTION	OUTRIGHT SALE PRICE	EXCHANGE SALE PRICE
43403	GRAPHIC ARTS & INDUSTRIAL	\$2,270.00	\$942.00

10 INCH RACKS

PART NO.	DESCRIPTION	OUTRIGHT SALE PRICE	EXCHANGE SALE PRICE
<u>DEVELOPER RACK</u>			
43840	ALL-PRO 2000M, HS, NDT;FILM QUICK CT	\$720.00	\$525.00
<u>FIXER RACK</u>			
43845	ALL-PRO 2000M, HS, NDT, FILM QUICK CT	\$665.00	\$525.00
<u>WASH/DRY RACK</u>			
43850	ALL-PRO 2000M, HS, NDT, FILM QUICK 10, A/T 2020	\$1,280.00	\$953.00

FILM TRANSPORT (8 INCH)

PART NO.	DESCRIPTION	OUTRIGHT SALE PRICE	EXCHANGE SALE PRICE
43855	DEVELOPER	\$685.00	\$525.00
43860	FIXER	\$620.00	\$453.00
43865	WASH/DRY	\$1,125.00	\$809.00

ACCENT CCD DIRECT DIGITAL INTRAORAL SENSOR

PART NO.	DESCRIPTION	PRICE
B5130	DISPOSABLE INFECTION CONTROL SHEATHS SIZE 1 (BOX OF 500)	\$75.00
B5230	DISPOSABLE INFECTION CONTROL SHEATHS SIZE 2 (BOX OF 500)	\$75.00
B5030	SENSOR STORAGE HANGER	\$10.00
B5030-10	SENSOR STORAGE HANGER (BOX OF 10)	\$100.00
B5040	CONTROL MODULE HOLSTER WITH BELT	\$20.00
B5040-10	CONTROL MODULE HOLSTER WITH BELT (BOX OF 10)	\$200.00
B5027	USB CONNECTION CABLE (5M LENGTH)	\$9.95
B3554	USB CONNECTION CABLE (2M LENGTH)	\$6.95
B5080	A/C BATTERY CHARGER	\$150.00
B5085	STANDARD BATTERY CHARGER (REPLACEMENT)	\$50.00
B5085-5	STANDARD BATTERY CHARGER (REPLACEMENT) (BOX OF 5)	\$250.00
B5075	ADAPTER BLUETOOTH TO USB	\$150.00
B5075-5	ADAPTER BLUETOOTH TO USB (BOX OF 5)	\$750.00
B5013	RECHARGEABLE BATTERIES (PACK OF 4)	\$20.00
B5013-10	RECHARGEABLE BATTERIES (10 PACKS OF 4)	\$200.00

AP SCANX

PART NO.	DESCRIPTION	PRICE
A8225	ARCH UPPER COVER ASSY FOR MODEL A8200	\$220.00
74362	ARCH UPPER COVER ASSY FOR MODELS B3810, B3810P, A8100	\$293.20
B2835	ARCH LOWER COVER ASSY FOR MODELS B2800P, B2800PV, B2800F, B2800FV	\$175.00
73767	ARCH LOWER COVER ASSY FOR MODELS A8200, B3810, B3810P, A8100	\$103.00
73248-0	BARRIER ENVELOPES, SIZE 0 (QTY. 100)	\$14.00
73248-1	BARRIER ENVELOPES, SIZE 1 (QTY. 100)	\$14.00
73248-2	BARRIER ENVELOPES, SIZE 2 (QTY. 300)	\$38.00
73248-3	BARRIER ENVELOPES, SIZE 3 (QTY. 100)	\$18.00
73248-4	BARRIER ENVELOPES, SIZE 4 (QTY. 50)	\$18.00
74190-8	CASSETTE 8"x10"	\$125.00
74190-10	CASSETTE 10"x12"	\$150.00
74190-14	CASSETTE 14"x17"	\$205.00
B2010	CLEANING SHEET, PACKAGE OF 12	\$60.00
B2020	CLEANING SHEET, PACKAGE OF 25	\$100.00
B2982F	DUST COVER; PACKAGED SCANX 14 IN-COUNTER	\$66.95
B2982	DUST COVER; PACKAGED SCANX 14 PORTABLE	\$66.95
73096	LINE CORD	\$23.45
74365	PCB; STEPPER MODEL CONTROL ASSY	\$90.65
73445-0	PHOSPHOR STORAGE PLATE (PSP), SIZE 0 (QTY. 2)	\$44.00
73445-1	PHOSPHOR STORAGE PLATE (PSP), SIZE 1 (QTY. 2)	\$44.00
73445-2	PHOSPHOR STORAGE PLATE (PSP), SIZE 2 (QTY. 4)	\$88.00
73445-3	PHOSPHOR STORAGE PLATE (PSP), SIZE 3 (QTY. 2)	\$44.00
73445-4	PHOSPHOR STORAGE PLATE (PSP), SIZE 4 (QTY. 2)	\$60.00
73578-57	PHOSPHOR STORAGE PLATE (PSP), TMJ 5"x7" (QTY. 1)	\$489.25
73578-6	PHOSPHOR STORAGE PLATE (PSP), 15CM X 30CM (QTY. 1)	\$795.00
73578-5	PHOSPHOR STORAGE PLATE (PSP), 5"x12" (QTY. 1)	\$775.00
73578-8	PHOSPHOR STORAGE PLATE (DENTAL), 8"x10" (QTY. 1)	\$850.00
73578-8M	PHOSPHOR STORAGE PLATE (MEDICAL), 8"x10" (QTY. 1)	\$850.00
73578-10	PHOSPHOR STORAGE PLATE 10"x12" (QTY. 1)	\$995.00
73578-14	PHOSPHOR STORAGE PLATE 14"x17" (QTY. 1)	\$1,250.00
73566-0	PLATE GUIDE, SIZE 0 (QTY. 1)	\$24.70
73566-1	PLATE GUIDE, SIZE 1 (QTY. 1)	\$24.70
73566-2	PLATE GUIDE, SIZE 2 (QTY. 1)	\$24.70
73566-3	PLATE GUIDE, SIZE 3 (QTY. 1)	\$24.70
B2160-8M	PLATE PROTECTOR 8"x10"	\$15.00
B2160-10	PLATE PROTECTOR 10"x12"	\$18.00
B2160-14	PLATE PROTECTOR 14"x17"	\$25.00
73470	PLATE TRANSFER BOX	\$25.00
73097	POWER SUPPLY	\$133.30
73471	POWER SUPPLY WITH LINE CORD	\$156.70
73050	TRANSPORT ARCH ASSY	\$3,500.00
B2750	TRANSPORT ARCH ASSY FOR MODELS B2100, B2100V	\$3,605.00

NOTES



Phone: 516-433-7676 / Fax: 516-433-7683
www.allproimaging.com