



ScanX[®] 12 Portable

The Only Truly Portable CR System. Just Right for the Equine Practice.

You can rely on ScanX 12 Portable for fast, detail-rich digital images on the spot – wherever you work. With ScanX Computed Radiography you are able to view high resolution digital images in just seconds on location; eliminating the need to ever return for retakes. Treatment can begin more quickly; benefitting both client and patient.

Equine Veterinarians choose the ScanX 12 Portable because it is the most durable, affordable and portable digital imaging system on the market.



ScanX 12 Portable

ScanX[®] 12 Portable

ScanX performs to the highest standards and gives you flexibility in the field all while eliminating the hassles of film-based imaging.

- Weighs only 37 lbs. and packed in a custom, durable soft case.
- ScanX utilizes phosphor storage plates (PSP) which replace film and are reusable thousands of times.
- Make a diagnosis on the spot; share with your clients anytime, anywhere.
- Easily transfer your existing radiographic technique to digital.
- Eliminate the need to ever return for retakes.



High-resolution digital images are available for viewing in seconds.

Specifications

Dimensions in Softcase:	17"W x 17"D x 16 ½"H
Weight:	37 lbs.
Electrical:	115 V, 60 Hz
Laser Classification:	Class 1 Laser Product. Compliance with FDA HHS 21 CFR 1040.10 and IEC 60825-1
Imaging Media:	Phosphor Storage Plates: 8" x 10", 10" x 12"
System Includes:	Power Supply and USB Cable
Software:	Sold separately

Computer Recommendations

CPU Speed:	2.8 GHz Pentium IV or higher.
Operating System:	Windows 2000 Professional with Service Pack 4 or later or Windows XP Professional with Service Pack 2
RAM:	2 GB
CD ROM Drive:	CD/RW/DVD
Hard Drive:	120 GB
Monitor:	17" color monitor (1024x768 or Higher)
Video Display Adapter:	256 MB RAM
Keyboard / Mouse:	Standard keyboard and mouse
USB Port:	USB 2.0



1295 Walt Whitman Road, Melville, New York 11747-3062
800-247-8324 • Fax: 516-433-7683 • www.allproimaging.com