



Provecta™ V X-Ray

For your dental radiographs, benefit from the high quality and versatile operation of the ProvectaV X-Ray system

The ProvectaV exhibits unmatched stability. Constructed of a high strength aluminum alloy, the X-Ray arm resists the innate tendency of the tubehead to drift, resulting in extremely precise positioning. ProvectaV is truly versatile. All exposure values are pre-programmed—simply select from any of the Veterinary specific icons. Settings may also be manually adjusted to suit your individual preferences. Custom settings can be stored in memory and exposure time can be set in seconds or pulses.

With ProvectaV you will obtain excellent results; digital or film: Whether you use film or require ultra-short digital exposure times, you can switch instantly between both mediums by pressing a single icon. No guesswork. All settings are automatically programmed.

Specifications

Electrical:	108 - 132V, 60 Hz
X-Ray Tube:	70 kVp, 8mA; 360° tubehead rotation
Main Circuit Breaker:	10 Amp/120V
Focal Distance:	8"/20 cm
Focal Spot:	8 mm (IEC 336)
Total filtration:	2 mm Aluminum Equivalent
Duty Cycle:	1:32
Exposure Range:	0.02 sec. / 1.0 pulses to 3.2 sec. / 192.0 pulses



ProvectaV Wall Mount



ProvectaV Mobile