



Specifications

- Dimensions in case: 15 1/2 W x 18 L x 14 H
- Weight: 44 lbs. (without battery) 46 lbs. (with battery)
- Electrical: 50/60 Hz | 110-240 Vac. 24 Volt DC power port
- Laser Classification: Class I Laser Product. Compliance with FDA HHS 21 CFR 1040.10 and IEC 60825-1.
- Imaging Media: Phosphor Storage Plates, any brand. All standard and custom sizes.
- System Includes: Power Supply, DC Power Adapter and USB Cable.
- Part Numbers: D5000K-D (ScanX Discover w/o battery)
D5000K-DB (Scanx Discover w/o battery)
D5000K-DR (ScanX Discover HR w/o battery)
D5000K-DRB (ScanX Discover HR w/ battery)