



Specifications

Dimensions in Softcase:	15 ½ W x 18 L x 14 H
Weight:	43 lbs. (without battery) 45 lbs. (with battery)
Electrical:	50/60 Hz 110-240 Vac, 24 volt DC power port
Laser Classification:	Class I Laser Product Compliance with FDA HHS 21 CFR 1040.10 and IEC 60825-1
Imaging Media:	Phosphor Storage Plates, all brands. All standard and custom sizes
System Includes:	Power Supply, DC Power Adapter and USB Cable
Software:	Compatible, third party software recommended (not included)
Part Numbers:	D5000-S (ScanX Scout w/o battery) D5000-SB (ScanX Scout w/ battery)

Computer Recommendations

CPU Speed:	3.0 GHz Intel CORE2
Operating System:	Microsoft Windows XP Professional (for a 32-bit Intel processor with Service Pack 2 or later or for a 64-bit Intel x64 processor) Microsoft Vista Business, Enterprise, or Ultimate (for a 32-bit Intel processor with Service Pack 1 or later or for a 64-bit Intel x64 processor with Service Pack 1 or later)
RAM:	2 GB
CD ROM Drive:	CD/RW/DVD
Hard Drive:	500 GB
Optical Drive:	Optical drive that is CD Rom compatible
Monitor:	SVGA 20", 1024x768 or higher resolution, contrast ratio 450:1, .22 dot pitch recommended
Video Display Adapter:	32 MB RAM
Keyboard / Mouse:	Standard keyboard and mouse
USB Port / Version:	USB 2.0 or later