



Specifications

Dimensions: 15 ½ W x 18 L x 14 H
Weight: 40 lbs. (without battery) 42 lbs. (with battery).
Electrical: 50/60 Hz 110-240 Vac, 24 volt DC power port.
Laser Classification: Class 1 Laser Product. Compliance with FDA HHS 21 CFR 1040.10 and IEC 60825-1.
Imaging Media: Phosphor Storage Plates, all brands. All standard and custom sizes.
System Includes: Power Supply, DC Power Adapter and USB Cable.
Software: Compatible, third party software recommended (not included).
Part Numbers: D5000K-T (ScanX Trek w/o battery)
D5000K-TB (ScanX Trek w/ battery)

Computer Recommendations

CPU Speed: 3.0 GHz Intel CORE2
Operating System: Microsoft Windows XP Professional (for a 32-bit Intel processor with Service Pack 2 or later or for a 64-bit Intel x64 processor) Microsoft Vista Business, Enterprise, or Ultimate (for a 32-bit Intel processor with Service Pack 1 or later or for a 64-bit Intel x64 processor with Service Pack 1 or later).
RAM: 2 GB
CD ROM Drive: CD/RW/DVD
Hard Drive: 500 GB
Optical Drive: Optical drive that is CD Rom compatible
Monitor: SVGA 20", 1024x768 or higher resolution, contrast ratio 450:1, .22 dot pitch recommended.
Video Display Adapter: 32 MB RAM
Keyboard / Mouse: Standard keyboard and mouse
USB Port/Version: USB 2.0 or later